

PERSONA<sup>®</sup>  
MEDICAL



Invisible never looked so good™

Ole



THE WORLD'S  
SMALLEST HEARING  
INSTRUMENT

with **AirTouch™** Memory

Invisible never looked so good!

## Open Fit, Back In The Ear Where It Belongs

From the inventor of Open Tube Technology comes Olé, the first open fit completely in canal hearing device with AirTouch™ Multi-Memory.

Olé is one of the finest hearing instruments ever developed. The Olé hearing instrument is the first open fit, completely in the canal hearing device with AirTouch™ Multi-Memory. Amazing sound quality, noise reduction and virtual invisibility make Olé exactly what you desire from a hearing device.

Olé products are also available in other traditional styles from custom fit to behind-the-ear instruments with additional features appropriate for various hearing loss.



Olé

Invisible never looked so good!



Olé

## The Perfect Fit

The human voice is one of the most complex sounds. The ability to understand speech, processing millions of computations per second, is the essence of digital technology. Bringing the sounds you need into your hearing range is perfectly attained through custom computerized programming by your hearing care professional.

Programming your Olé allows for the most accurate fitting for your hearing loss and unprecedented enhancements.

Invisible never looked so good!


- Digital / Programmable
- No Feedback / No Whistle
- Tinnitus Masker Option
- No Tubes, Tips or Molds
- Retains Natural Ear Resonance
- AirTouch™ Multi Memory
- Easy Placement
- High Frequency Power
- No Echo
- Secure For Active Patients
- Better Comfort / No Tickle
- Retains Natural Directionality
- Less Perspiration Damage
- Minimize Maintenance
- Wax Barrier

Invisible never looked so good!

## AirTouch Multi-Memory

AirTouch Memory allows you to select between memories with the greatest of ease. Simply tap the outer part of your ear canal and the sound pressure triggers AirTouch. The Olé does not require buttons or switches and does not need to be touched to change memories.

With your fingers together, lightly tap the ear canal area. This will toggle successively from the first memory through the fourth, and then back to the first.



Invisible never looked so good!

## Invisible Never Looked so Good!

Olé fits comfortably inside your ear canal and is virtually invisible allowing you to feel comfortable and confident. Don't be fooled by its small size however; the Olé is loaded with the latest hearing aid features such as our state of the art Noise reduction system, the most advanced feedback canceller and AirTouch Multi-Memory.

Consult your hearing care professional today and start experiencing Olé for yourself.



Invisible never looked so good!



Olé

800.789.6543  
170 N. Cypress Way • Casselberry, FL 32707  
Info@PersonaMedical.com  
PersonaMedical.com



BR / OLE/ 3 © 2011 Persona Medical

P E R S O N A<sup>®</sup>  
M E D I C A L